

How to clean a CGM transmitter



First, make sure you have your watertight tester. Properly connecting the tester to the transmitter ensures that fluids do not come in contact with the connector pins inside the transmitter.

Before you prepare to clean your transmitter, ensure you have the tester available.



Gather these supplies

- Mild liquid soap
- Soft-bristled small toothbrush
- Cotton swab
- Container
- Clean lint-free dry cloths



If you will also be removing adhesive residue from your transmitter, you also need a medical adhesive remover such as Zoff®, Unisolve®, Remove®, and cotton swabs. We recommend that you do this after you've removed the transmitter, and before your next sensor insertion.

Cleaning your CGM transmitter

1. Wash your hands thoroughly.
2. Attach the tester to the transmitter by lining up the flat side of the tester with the flat side of the transmitter. Push the tester into the transmitter until you hear a click. When properly connected, the green light flashes.
3. Holding the tester, use the remover solution wipe to clean the transmitter.
4. Rinse the transmitter under room temperature tap water for at least one minute.
5. Prepare a mild liquid soap solution using 1/2 teaspoon of mild liquid soap per litre of room temperature tap water. With the tester still attached, submerge in the mild liquid soap
6. Dry the transmitter still connected to the tester with a clean, dry cloth.
7. Disconnect from the transmitter by gently squeezing the arms of the tester and pulling away from the transmitter.



Indicated sensor insertion locations vary by age. See user guide.