


MAKE THE CHANGE

For better control and to maintain skin health, change your infusion set and site every 2-3 days.




PROBLEMS THAT CAN OCCUR AFTER 3 DAYS

Why it is important to change your infusion set every 2-3 days.

Poor Insulin Absorption

Longer wear may result in poor insulin absorption and lead to higher glucose levels.^{1,2}




In addition, insulin loses its strength over time³, especially with heat and agitation. As this happens, insulin becomes less effective, and controlling your blood glucose becomes more difficult. Therefore, it is critical to change your insulin reservoir every 2-3 days as well.

Site irritation

People who use an infusion set for longer than three days are four times more likely to experience site irritations¹.

■ Itching ■ Bruising ■ Swelling ■ Pain




More set problems, such as poor adhesion and leakage start to occur the longer the set is used¹.

ROTATE INFUSION SITES REGULARLY


Each time you change your infusion set you should also rotate the site.

Re-using the same site can lead to the fatty tissue under your skin becoming hard and bumpy (lipohypertrophy). Insulin is not absorbed evenly in these areas and can lead to unexpected high or low blood glucose levels.

Below are two examples of patterns use to help you rotate your sites.



Imagine a clock drawn on your abdomen. Rotate infusion sites by starting at the 12 o'clock position and use the site at 3 o'clock, then 6 o'clock, and so on.



Imagine an 'M' or 'W' pattern approximately 10 cm long on either side of your belly button. Start at the end of one letter and proceed through the letter, rotating the infusion set at each corner.

It does not matter what rotation method (or combination of methods) you use as long as you rotate the infusion sites and avoid the 5 cm area surrounding your belly button.

CHOOSING THE RIGHT INFUSION SET

Your choice of infusion set should be based on your lifestyle, age, build and clinical recommendations.

Plastic Cannula

MINIMED® QUICK-SET™

COMFORT AND EASE

Our most popular infusion set combining ease of use with maximum comfort. Its virtually painless insertion with the Quick-Serter™ device, even in hard-to-reach areas, makes it a great choice for children, adults of average to large build and those who have limited dexterity or are less active.

FEATURES:

- 90° soft cannula, consistent insertion depth for optimal insulin absorption
- Comfortable, virtually painless insertions with the Quick-Serter™
- Convenient disconnection at the infusion site



MINIMED® MIO™

ALL-IN-ONE CONVENIENCE

The newest MiniMed® all-in-one infusion set, integrating the set and the insertion device in one single unit that comes in a variety of colours. Easy to insert and carry, MiniMed® mio™ is a great choice if you're looking for a portable infusion set, are less active and looking for alternative site usage.

FEATURES:

- Integrated inserter, for virtually painless insertion
- 90° soft Teflon cannula
- Available in transparent, blue or pink



CHOOSING THE RIGHT INFUSION SET

Your choice of infusion set should be based on your lifestyle, age, build and clinical recommendations.

MINIMED® SILHOUETTE™

COMFORT AND FLEXIBILITY

The MiniMed® infusion set offering the most flexibility, with a variable insertion angle from 20° to 45° enabling you to adjust the depth of insertion. Combined with Silhouette's thin profile it makes it particularly comfortable for lean children or adults, very active people or pregnant women (from 2nd trimester).

FEATURES:

- 90° soft Teflon cannula
- Adjustable depth of insertion, from 20° to 45° angle
- Virtually painless insertion with Sil-Serter.



Steel Cannula

MINIMED® SURE-T™

ENHANCED SECURITY AND SIMPLICITY

Innovatively designed to provide enhanced security for continuous insulin delivery. Its kink-free steel needle, the thinnest available in the MiniMed® portfolio, combined with an extra disconnection site 10 cm away from the infusion site provides extra security against clogging and dislodging. It is particularly suitable for infants and children, pregnant women (up to 2nd trimester) or people allergic to Teflon®.



FEATURES:

- Ultra-thin steel needle
- Easy 90° manual insertion
- Disconnection 10 cm away from the infusion site






MINIMED INFUSION SETS

PRODUCT CODES

			LENGTH		MEDTRONIC CODE
NDSS Code	PRODUCT NAME		Cannula/ Needle (mm)	Tubing (cm)	MiniMed® connector (for Paradigm 511, 512/712, 515/715, 522/722, Veo & MiniMed series pumps)
564	MiniMed® Quick-set™ (10 per box) 		6	60	MMT-399A
659			6	80	MMT-387A
563			6	110	MMT-398A
562			9	60	MMT-397A
660			9	80	MMT-386A
561			9	110	MMT-396A
640	MiniMed® Silhouette™ (10 per box) 		13	46	MMT-368A
588			13	60	MMT-381A
641			13	80	MMT-383A
589			13	110	MMT-382A
569			17	60	MMT-378A
642			17	80	MMT-384A
568			17	110	MMT-377A
647	MiniMed® Sure-T® (10 per box) 		6	60	MMT-864A
650			6	80	MMT-866A
648			8	60	MMT-874A
651			8	80	MMT-876A
649			10	60	MMT-884A
652			10	80	MMT-886A
643	MiniMed® Mio® (10 per box)		6	46	MMT-921A
644		Pink	6	60	MMT-923A
645			6	80	MMT-925A
643.1			6	46	MMT-941A
644.1		Blue	6	60	MMT-943A
645.1			6	80	MMT-945A
645.2		Clear	6	80	MMT-965A
646			9	80	MMT-975A
705	MiniMed® insulin reservoir (10 per box)		1.8 mL	-	MMT-326A
706			3.0 mL	-	MMT-332A
n/a	MiniMed® Quick-serter™ insertion device (purchase at www.medtronic-diabetes.com.au , or Medtronic Customer Service on 1800 777 808, option 1)				MMT-395QS
n/a	MiniMed® Sil-serter™ insertion device (purchase at www.medtronic-diabetes.com.au , or Medtronic Customer Service on 1800 777 808, option 1)				MMT-385

MINIMED INFUSION SETS AVAILABLE IN LUER LOCK

NDSS Code	PRODUCT NAME	LENGTH		MEDTRONIC CODE
		Cannula/ Needle (mm)	Tubing (cm)	Luer-lock connector
579	MiniMed® Quick-set™ (10 per box) 	6	60	MMT-393
577		9	60	MMT-392
640	MiniMed® Silhouette™ (10 per box) 	13	46	MMT-368
588		13	60	MMT-381
653	MiniMed® Sure-T® (10 per box) 	6	60	MMT-863
654		8	60	MMT-873
658		10	80	MMT-885

Infusion Set Fill Amounts

MIO®

Fill Amount

6mm – 0.300 u

9mm – 0.500 u

QUICK-SET™

Fill Amount

6mm – 0.300 u

9mm – 0.500 u

SILHOUETTE™

Fill Amount

13mm – 0.700 u

17mm – 0.700 u

SURE-T®

When using Sure-T
you are not required
to fill the cannula.

Infusion sets are available via Pharmacy and for non-Eligible NDSS please visit Medtronic Diabetes eShop:

<https://eshop.medtronic-diabetes.com.au/>
or contact

Medtronic Customer Service 1800 777 808 (Option 1)

* Visit <http://www.ndss.com.au> to see the NDSS eligibility criteria.

TIPS FOR PREVENTING INFECTIONS

- Always wash your hands before handling infusion sets and sites
- Thoroughly prepare the infusion site (as instructed by your healthcare provider)
- Avoid touching the needle, cannula or other parts that come in to contact with insulin
- Check the site frequently for any problems
- Change infusion sets and sites every 2-3 days
- Do not keep sets if the packaging seal is broken or the expiration date has passed

POST-REMOVAL SKIN CARE

Adhesive-removal products can be used to help get rid of any left behind adhesive. You can also apply an antiseptic cream to help prevent infections or a moisturiser to prevent any dryness.

DON'T FORGET THE RESERVOIR

Insulin reservoirs are important parts of your insulin pump system. It is also important to change these reservoirs every 2-3 days, particularly if they have been exposed to excessive heat or sunlight.⁴

ALWAYS READ THE INSTRUCTIONS FOR USE.

References: 1. Schmid V, Hohberg C, Borchert M, Forst T, Pfützner A. Pilot study for assessment of optimal frequency for changing catheters in insulin pump therapy-trouble starts on day 3. *JDiabetesSciTechnol*. 2010Jul1;4(4):976-82. 2. ThethiTK, RaoA, KawjiH *et al*. Consequences of delayed pump infusion line change in patients with type 1 diabetes mellitus treated with continuous subcutaneous insulin infusion. *J Diabetes Complications*. 2010 Mar-Apr;24(2):73-8. 3. NovoLog®. Princeton, NJ: Novo Nordisk Inc.; 2007 4. Novo Nordisk Pharmaceuticals. NovoRapid-NovoMix Product Information. Baulkham Hills NSW. July 2014.

Medtronic

Medtronic Australasia Pty Ltd
2 Alma Road,
Macquarie Park, NSW 2113
Australia
Tel: +61 2 9857 9000
Fax: +61 2 9889 5167
Toll Free: 1800 668 670

E-Mail: australia.diabetes@medtronic.com
Facebook: www.facebook.com/MedtronicDiabetesAUS
Twitter: @DiabetesANZ
YouTube: Medtronic Diabetes ANZ

APPROVAL #7773-052020 © Medtronic 2020
All Rights Reserved. Printed in Australia 05/2020

medtronic-diabetes.com.au