

Assessment & Progress Report Guide

Medtronic CareLink™ Personal software



Representative patient profile

- These are your selected data date ranges. Data A (blue) is your most recent selected date range i.e. the past 2 weeks. Data B (orange) is your previous data range to make comparison easy.
- This graph represents all of your sensor readings over 24hrs (midnight to midnight). The darker solid areas represent where most of your SG (sensor glucose) reading fell during this period. The coloured dotted line represents your high and low variability in this time period. A narrower band would indicate less variability.
- The dark dotted line represents your average SG levels for data A only.
- This is your time in range (TIR) data between 3.9 - 10 mmol/L. You can easily compare this information between columns A and B to easily track any changes.
- If your MiniMed™ 780G system exits SmartGuard™ then the reasons for exit will be listed here. On your sensor trace you will see these exits as numbered black boxes.
- These figures show how long you spent in the SmartGuard™ feature. You can also review your number of high and low alarms. Your glucose management indicator (GMI) may be similar to your laboratory HbA1c but it is an approximate so could be slightly different.
- This is your daily average blood glucose and calibration data.
- This is your insulin use data. How much you are using altogether, how much is basal and how much is bolus. You can also check how much Auto Correction insulin you are getting and how often you have been changing your infusion set and reservoir.
- This section covers your meals. You can review how many meals and carbohydrates you are eating on average per day.

Personalised diabetes insights anytime, anywhere

Medtronic

Medtronic Australasia Pty Ltd
2 Alma Road
Macquarie Park, NSW 2113, Australia

Tel: +61 2 9857 9000
Fax: +61 2 9889 5167
Toll Free: 1800 777 808

E-Mail: australia.diabetes@medtronic.com
Facebook: www.facebook.com/MedtronicDiabetesAUS
YouTube: [MedtronicDiabetesANZ](https://www.youtube.com/MedtronicDiabetesANZ)

<https://www.medtronic-diabetes.com.au/>

Always read the instructions for use.

For detailed information regarding the instructions for use, indications, contraindications, warnings, precautions, and potential adverse events, please consult the device manual. For further information, contact your local Medtronic representative. © 2021 Medtronic Australasia Pty Ltd. All Rights Reserved. Please note: In contacting the Diabetes Toll Free, your personal and health information may be disclosed to an operator located outside Australia. 10312-102021