

InPen™ app Compatibility Guide

Medtronic

InPen™ App version 7.0.0

Compatible OS versions

iOS	14, 15, 16, 17
------------	----------------

Android	11, 12, 13
----------------	------------

Incompatible devices

Device	Model Number
---------------	---------------------

All tablets	
-------------	--

LG Stylo Series	
-----------------	--

Additional notes

Bluetooth® performance varies due to differences in phone hardware and operating systems. The InPen App is not recommended on Jailbroken or Rooted devices.

For InPen™ app, all devices are either fully compatible or unverified unless explicitly listed in the table of “Incompatible devices” above.