



# 25 REASONS FOR **UPGRADE**

# YOU'VE BEEN USING THE SAME INSULIN PUMP FOR NEARLY FOUR YEARS SEE 25 REASONS TO UPGRADE!



## ADVANCED PROTECTION

**1. SMARTGUARD™** mimics some of the functions of a healthy pancreas by automatically stopping insulin delivery when sensor glucose is predicted to approach a low and resuming delivery when levels recover.

**2. JUMP IN THE POOL WITH CONFIDENCE** The MiniMed® 640G is waterproof to 3.6 metres for up to 24 hours, so you don't have to worry when it comes to bathing or getting in the water.

**3. EXPERIENCE THE PRECISION OF BOLUS WIZARD™** Bolus Wizard™ calculates insulin doses automatically so you can manage your insulin more precisely.

**4. GET CONTINUOUS PEACE OF MIND WITH CGM** Continuous Glucose Monitoring tracks your glucose levels 24 hours a day to provide a complete picture of your glucose levels and give you early warnings.

## GREATER CONVENIENCE

**5. SEE EASIER WITH A NEW COLOUR SCREEN** Experience the clearer vision of a colour screen with improved resolution and 5 levels of brightness variability

**6. ACTIVE INSULIN AT A GLANCE** When used with CGM, the MiniMed 640G displays your blood glucose and active insulin levels providing a full picture.

**7. PRE-SET BOLUS VALUES** Where once there were no pre-set bolus values with your Veo the MiniMed640G can hold up to 8 values giving you greater convenience at meal time with pre-set functions for breakfast, lunch, dinner, snack-time and four other boluses of your choice.

**8. ADJUSTABLE BOLUS SPEED** The MiniMed 640G allows for a choice of two bolus speeds.

**9. LINKABLE DEVICES** The advantage of having

an additional 3 linkable devices allows you to connect multiple Bayer meters. Giving you the convenience of leaving a meter at work, in the car, gym bag and at home without having to worry about forgetting it.

## 10. INFUSION SET-CHANGE REMINDERS

With many users having a busy schedule it is often easy to forget how many days it has been since your last set change – the MiniMed 640G provides a set change reminder.

**11. BOLUSES FROM METER** You can achieve greater convenience with the ability to Bolus from the meter without having to touch your pump – allowing you to bolus discretely while out at dinner with friends.

**12. ALARMS** No longer do you have to choose between audio or vibration alerts, you can now have both.

**13. VOLUME** The MiniMed 640G allows you to adjust the volume of your alerts. You can now set it to a level twice as loud as the Veo

**14. ACTIVE INSULIN DISPLAYED ON HOME SCREEN** You can now see the active amount of insulin you have on board by glancing at your pump

**15. UPLOADING DATA TO CARELINK** While uploading your data into CareLink the MiniMed 640G delivers insulin normally – Veo normally suspends insulin delivery until upload is complete.

**16. NEW METER-CONTOUR NEXT LINK** The CONTOUR® NEXT LINK 2.4 has a modern compact design with quick and easy readings after a five second countdown. The added Second-Chance™ sampling tells you to apply more blood if your first sample was not sufficient.

## IMPROVED DESIGN

**17. WATERPROOF** The new MiniMed 640G is waterproof for up to 3.6m for 24 hours†, giving you the freedom to enjoy a swim in the warmer

months of the year – the Veo model pump has a 'water-tight' seal of 1 meter for up to 30 minutes

**18. DISPLAY QUALITY** Providing a full colour display of 65,000 colours – the MiniMed 640G provides a greater user experience from Veo's Monochrome display.

**19. DISPLAY RESOLUTION** The MiniMed 640G has more than double the amount of pixels – 320 X 240 compared to the Veo with a resolution of 120 X 32 pixels.

**20. SCREEN DIMENSIONS** With an increased area of around 4cm<sup>2</sup> the MiniMed 640G allows for an easier to read display with greater ease of use.

**21. BUTTON BACKLIGHT** New to Medtronic pump technology the MiniMed640G has an automatic button backlight function which provides you with the ability to interact easily with your pump in the dark.

**22. ADJUSTABLE BACKLIGHT** The 30 second backlight previously available on the VEO has been upgraded to an LED display which is adjustable from 30 seconds – 2 minutes

**23. BRIGHTNESS AUTO ADJUST** The MiniMed 640G's brightness auto adjust allows the user to easily interact with the pump when outdoors. Where there was once no brightness variability with the VEO the 640G also provides a manually adjustable scale of 1-5.

**24. DEVICE ORIENTATION** Through changing the orientation of the pump from landscape to portrait makes the product more ergonomic by allowing both the left and right hand dominant to use the pump with ease.

**25. NAVIGATION** Providing 5 navigation buttons with the MiniMed 640G allows the user to navigate through a variety of options easily – four directional buttons and the menu button allow for a quick return to the menu page.

**Medtronic**