

# InPen™ App

## Compatibility Guide

Medtronic

### InPen™ App version 7.2.0

#### Compatible OS versions

<b>iOS</b>	14, 15, 16, 17, 18
<b>Android</b>	11, 12, 13, 14

#### Incompatible devices

<b>Device</b>	<b>Model Number</b>
All tablets	
LG Stylo Series	
Huawei Brand	All – Blocklisting in MAY 2025

#### Additional notes

Bluetooth® performance varies due to differences in phone hardware and operating systems. The InPen App is not recommended on Jailbroken or Rooted devices.

*For InPen™ app, all devices are either fully compatible or unverified unless explicitly listed in the table of "Incompatible devices" above.*