

# WHICH INFUSION SET IS RIGHT FOR YOU?


## STEEL

- Reactions to plastic cannula
- History of bent cannula
- Pregnant (up to 2nd trimester)
- Small child

**90 degree insertion**

### SURE-T® INFUSION SET

Simple steel needle infusion set




## PLASTIC CANNULA

- Less active
- Alternative site usage
- Limited dexterity

**90 degree insertion**

### QUICK-SET™ INFUSION SET

Our most popular infusion set




- Less active
- Alternative site usage

**90 degree insertion**

### MIO® INFUSION SET

Convenient "all in one" infusion set



- Active
- Lean
- Muscular
- Pregnant (from 2nd trimester)

**30-45 degree insertion**

### SILHOUETTE™ INFUSION SET

