

Medtronic Infusion Sets



Medtronic Extended™

Long Wear



MiniMed™ Sure-T™

Secure



MiniMed™ Quick-Set™

Comfortable



MiniMed™ Silhouette™

Adjustable Insertion Angle

User Needs & Special Circumstances

Suited to all body types. Fewer set changes¹ Dexterity issues and needle anxiety²

Soft cannula allergies
History of bent cannula

Suitable for all body types.

Suitable for lean and active body type

Insertion Device

(Integrated)serter

Manual

MiniMed™ Quick-serter

Manual

Cannula/ Needle Type

Soft

Steel

Soft

Soft

Insertion Angle

90°

90°

90°

20-45°

Cannula Length

6mm

9mm

6mm

8mm

6mm

9mm

13mm

17mm

Cannula Fill

0.6 Units

0.6m Units

No cannula fill required

0.3 Units

0.5 Units

0.7 Units

0.7 Units

Tubing Length

60cm, 80cm

60cm, 80cm

60cm, 110cm

60cm, 110cm

Cannula Thickness (Gauge}

25G

29G

25G

25G

Change Frequency

Up to 7 days

1-2 days

2-3 days

2-3 days

* Vs the 3-day standard infusion set

1. Ilany J, et al. Abstract 416 - Clinical study of a new extended wear infusion set design. 13th ATTD International Conference. February 2020. (2) Medtronic Data on File. MiniMed™ Mio™ Advance Usability Study 10619721 DOC Version Patients interviewed in the US. 2020

2. Medtronic Data on File MiniMed™ Mio™ Advance Usability Study 10619721 DOC Version 8, 30 Patients interviewed in the US. 2020.

Important Information ALWAYS FOLLOW THE DIRECTIONS FOR USE.

This material does not replace or supersede the instructions for use. It should not be considered the exclusive source of information, and should be used in conjunction with the User Guide. See the User Guide for detailed information on the instructions for use, indications, contraindications, warnings, precautions, and potential adverse events. For further information, contact your local Medtronic representative and/or consult the Medtronic ANZ website at medtronic-diabetes.com.au

Medtronic

Medtronic Extended Wear™ Sets and Reservoir

Medtronic Code	NDSS code	Description	Cannula (mm)	Tubing (cm)
MMT-431AK	661	Medtronic Extended™ Infusion Set	6mm	60cm
MMT-432AK	662	Medtronic Extended™ Infusion Set	6mm	80cm
MMT-441AK	663	Medtronic Extended™ Infusion Set	9mm	60cm
MMT-442AK	664	Medtronic Extended™ Infusion Set	9mm	80cm
MMT-342G	707	Medtronic Extended™ Reservoir 3.0ml	3.0ml	



Medtronic MiniMed™ Infusion Set and Reservoir

Medtronic Code	NDSS code	Description	Cannula (mm)	Tubing (cm)
MMT864A	647	Medtronic™ MiniMed Sure-T	6mm	60cm
MMT-866A	650	Medtronic™ MiniMed Sure-T	6mm	80cm
MMT-874A	648	Medtronic™ MiniMed Sure-T	8mm	60cm
MMT-876A	651	Medtronic™ MiniMed Sure-T	8mm	80cm
MMT-399A	564	Medtronic™ MiniMed Quickset	6mm	60cm
MMT-397A	562	Medtronic™ MiniMed Quickset	9mm	60cm
MMT-398A	563	Medtronic™ MiniMed Quickset	6mm	110cm
MMT-396A	561	Medtronic™ MiniMed Quickset	9mm	110cm
MMT-381A	588	Medtronic™ MiniMed Silhouette	13mm	60cm
MMT-378A	569	Medtronic™ MiniMed Silhouette	17mm	60cm
MMT-377A	568	Medtronic™ MiniMed Silhouette	17mm	110cm
MMT-332A	706	Medtronic™ MiniMed Reservoir	3.0ml	
MMT-395QS		MiniMed™ Quick-serter™ insertion device (purchase at www.medtronic-diabetes.com.au , or Medtronic Customer Service on 1800 777 808, option 1)		

