










MiniMed™ 780G System

Simplified Menu Map



<p> HISTORY & GRAPH</p> <p>History Sensor Glucose Review Graph Time in Range</p>	<p> SMARTGUARD</p> <p>SmartGuard Checklist Temp Target SmartGuard Settings SmartGuard On and Off</p>	<p> SOUND & VIBRATION</p> <p>Silence Sensor Alerts Volume Sound Vibration Alert Settings Shortcut</p>
<p> RESERVOIR & SET</p> <p>New Reservoir Fill Cannula</p>	<p> INSULIN</p> <p>Bolus Basal Suspend/Resume Basal Delivery Delivery Settings Shortcut</p>	<p> BLOOD GLUCOSE</p> <p>BG</p>
<p> STATUS</p> <p>Suspend All Delivery SmartGuard Checklist Pump Sensor</p>	<p> PAIRED DEVICES</p> <p>Pair New Device Pair CareLink Sensor</p>	<p> SETTINGS</p> <p>Alert Settings Delivery Settings Device Settings</p>