










<div></div> <div>HISTORY & GRAPH</div> <div>History</div> <div>Sensor Glucose Review</div> <div>Graph</div> <div>Time in Range</div>	<div></div> <div>SMARTGUARD</div> <div>SmartGuard Checklist</div> <div>Temp Target</div> <div>SmartGuard Settings</div> <div>SmartGuard On and</div>	<div></div> <div>SOUND & VIBRATION</div> <div>Silence Sensor Alerts</div> <div>Volume</div> <div>Sound</div> <div>Vibration</div> <div>Alert Settings shortcut</div>
<div></div> <div>RESERVOIR & SET</div> <div>New Reservoir & Set</div> <div>New Reservoir Only</div> <div>New Set Only</div> <div>Fill Cannula</div>	<div></div> <div>INSULIN</div> <div>Bolus</div> <div>Basal</div> <div>Suspend All Delivery</div> <div>Delivery Settings shortcut</div>	<div></div> <div>BLOOD GLUCOSE</div> <div>BG</div>
<div></div> <div>STATUS</div> <div>Suspend All Delivery</div> <div>SmartGuard Checklist</div> <div>Pump</div> <div>Sensor</div>	<div></div> <div>PAIRED DEVICES</div> <div>Pair New Device</div> <div>Pair CareLink</div>	<div></div> <div>SETTINGS</div> <div>Alert Settings</div> <div>Delivery Settings</div> <div>Device Settings</div>

Medtronic Australasia Pty Ltd
Level 8, 11 Khartoum Road
Macquarie Park, NSW 2113, Australia

Tel: +61 2 9857 9000
Fax: +61 2 9889 5167
Toll Free: 1800 777 808

[Always read the instructions for use](#)

For detailed information regarding the instructions for use, indications, contraindications, warnings, precautions, and potential adverse events, please consult the device manual. For further information, contact your local Medtronic representative.
© 2025 Medtronic Australasia Pty Ltd. All Rights Reserved. Please note: In contacting the Diabetes Toll Free, your personal and health information may be disclosed to an operator located outside Australia.